

Data Sheet

Mixing / Magnetic Stirrers



C-MAG HS 7 digital IKAMAG®

New magnetic stirrer with heating and ceramic heating plate which offers excellent chemical resistance.

- Powerful motor for stirring quantities of up to 10 l (H₂O)
- Simultaneous digital display of target and actual temperatures via LCD display
- Directly connection for a PT 1000 temperature sensor enable a precise temperature control (included in delivery)
- Control accuracy in medium +/- 0,5 K (in combination with PT 1000)
- Display of actual temperature value in medium with 0,1 K resolution when using PT 1000 temperature sensor
- 3 operating modes to choose (standard, safe, adjustment protection)
- Fixed safety circuit of 550 °C
- Hot Top indicator >> hot surface warning to prevent burns!
- Digital error code display
- Elevated control panel for protection against leaking liquids

Accessories: H 38 Holding rod, H 44 Boss head clamp, RS 1 Set of magnetic stirring bars, IKAFLON® 50 Set (5 Pcs) slide round PTFE-coated Magnetic stirring bars, IKAFLON® 9 Set (5 Pcs) crown PTFE-coated Magnetic stirring bars, IKAFLON® 125 (1 Pc) beaker PTFE-coated Magnetic stirring bars, IKAFLON® 10 Set (5 Pcs) round PTFE-coated Magnetic stirring bars, IKAFLON® 15 Set (5 Pcs) round PTFE-coated Magnetic stirring bars, IKAFLON® 20 Set (5 Pcs) round PTFE-coated Magnetic stirring bars, IKAFLON® 25 Set (5 Pcs) round PTFE-coated Magnetic stirring bars, IKAFLON® 30 Set (5 Pcs) round PTFE-coated Magnetic stirring bars, IKAFLON® 40 Set (5 Pcs) round PTFE-coated Magnetic stirring bars, IKAFLON® 50 Set (5

Technical Data	
Number of stirring positions	1
Stirring quantity max. per stirring position (H ₂ O) [l]	10
Stirring quantity max. (H ₂ O) [l]	10
Motor rating input [W]	15
Motor rating output [W]	1.5
Speed display	scale
Speed range [rpm]	100 - 1500
Stirring bar length max. [mm]	80
Heat output [W]	1000
Heating rate (1 l H ₂ O in H15) [K/min]	5
Heating temperature range [°C]	50 - 500
Heat control	stepless
Heat control accuracy [±K]	1
Speed control	scale 1-10
Fixed safety circuit [°C]	550
Connection for ext. temperature sensor	PT1000
Control accuracy with sensor [±K]	0.5
Temperature constancy in medium [±K]	0.5
Accuracy of temperature measurement [K]	±0.2 + tolerance PT1000 (DIN IEC 751 Class A)
Temperature measurement resolution [K]	0.1
Set-up plate material	ceramic
Set-up plate dimensions [mm]	180 x 180
Dimensions (W x H x D) [mm]	215 x 105 x 330
Weight [kg]	5
Permissible ambient temperature [°C]	5 - 40
Permissible relative moisture [%]	80
Protection class according to DIN EN 60529	IP 21
Voltage [V]	230 / 120 / 100
Frequency [Hz]	50/60
Power input [W]	1020
Ident. No.	0003487001

Pcs) round PTFE-coated Magnetic stirring bars, IKAFLON® 80
 Set (5 Pcs) round PTFE-coated Magnetic stirring bars,
 IKAFLON® glass 25 Set (5 Pcs) round Magnetic stirring bars,
 IKAFLON® glass 30 Set (5 Pcs) round Magnetic stirring bars,
 IKAFLON® glass 40 Set (5 Pcs) round Magnetic stirring bars,
 IKAFLON® glass 50 Set (5 Pcs) round Magnetic stirring bars,
 IKAFLON® 15 Set (5 Pcs) power SmSo PTFE-coated Magnetic
 stirring bars, IKAFLON® 30 Set (5 Pcs) power SmSo
 PTFE-coated Magnetic stirring bars, IKAFLON® 50 Set (5 Pcs)
 power SmSo PTFE-coated Magnetic stirring bars, IKAFLON® 25
 Set (5 Pcs) slide round PTFE-coated Magnetic stirring bars,
 IKAFLON® 30 Set (5 Pcs) slide round PTFE-coated Magnetic
 stirring bars, IKAFLON® 40 Set (5 Pcs) slide round
 PTFE-coated Magnetic stirring bars, RS 2 Set of magnetic
 stirring bars, IKAFLON® 20 Set (5 Pcs) ellipse PTFE-coated
 Magnetic stirring bars, IKAFLON® 25 Set (5 Pcs) ellipse
 PTFE-coated Magnetic stirring bars, IKAFLON® 30 Set (5 Pcs)
 ellipse PTFE-coated Magnetic stirring bars, IKAFLON® 35 Set
 (5 Pcs) ellipse PTFE-coated Magnetic stirring bars, IKAFLON®
 40 Set (5 Pcs) ellipse PTFE-coated Magnetic stirring bars,
 IKAFLON® 50 Set (5 Pcs) ellipse PTFE-coated Magnetic stirring
 bars, IKAFLON® 60 Set (5 Pcs) ellipse PTFE-coated Magnetic
 stirring bars, IKAFLON® 70 Set (5 Pcs) ellipse PTFE-coated
 Magnetic stirring bars, IKAFLON® 10 Set (5 Pcs) cross
 PTFE-coated Magnetic stirring bars, IKAFLON® 20 Set (5 Pcs)
 cross PTFE-coated Magnetic stirring bars, IKAFLON® 25 Set (5
 Pcs) cross PTFE-coated Magnetic stirring bars, IKAFLON® 30
 Set (5 Pcs) cross PTFE-coated Magnetic stirring bars,
 IKAFLON® 38 Set (5 Pcs) cross PTFE-coated Magnetic stirring
 bars, IKAFLON® 37 Set (5 Pcs) bone PTFE-coated Magnetic
 stirring bars, IKAFLON® 54 Set (5 Pcs) bone PTFE-coated
 Magnetic stirring bars, TRIKA® 25 Set (5 Pcs) PTFE-coated
 Magnetic stirring bars, TRIKA® 40 Set (5 Pcs) PTFE-coated
 Magnetic stirring bars, TRIKA® 55 Set (5 Pcs) PTFE-coated
 Magnetic stirring bars, TRIKA® 80 Set (5 Pcs) PTFE-coated
 Magnetic stirring bars, IKAFLON® 67 (1 Pc) beaker
 PTFE-coated Magnetic stirring bars, IKAFLON® 74 (1 Pc)
 beaker PTFE-coated Magnetic stirring bars, IKAFLON® 103 (1
 Pc) beaker PTFE-coated Magnetic stirring bars, RSE Stirring bar
 remover, H 16 V Support rod, LABLIFT m



Quasar Instruments, LLC
4835 Centennial Blvd.
Colorado Springs, CO 80919

T. 877-535-8001 F. 719-574-2404 E. Sales@qi-us.com