

Ripette® pro

THE HANDY ONE

patent pending

Flexible: With the Ripette pro you can work flexibly and practically with 12 tip sizes, 10 volume adjustments and 120 program steps.

Clear: Volume adjustments and dosing volumes for different tips are shown in a vision panel on the Ripette pro.

Handy: A perfect ergonomic design makes an easy and comfortable use possible. With the easy-moving dosing button in combination with the low weight, long dosing series can be dispensed conveniently.

Compatible: The Ripette pro is compatible not only with PD-tips from Brand but also with Combitips and Combitips plus from Eppendorf and has been precisely aligned to Ritips professional from Ritter.

Independent: In addition to the handy design, the Ripette pro has been constructed consistent mechanical. The Ripette pro is maintenance-free and does not need any additional resources such as power supply or batteries.

		Program number									
Tip size	1	2	3	4	5	6	7	8	9	10	
0,05 ml	1	1,5	2	2,5	3	3,5	4	4,5	5	5,5	
0,10 ml	2	3	4	5	6	7	8	9	10	11	
0,20 ml	4	6	8	10	12	14	16	18	20	22	
0,50 ml	10	15	20	25	30	35	40	45	50	55	
1,00 ml	20	30	40	50	60	70	80	90	100	110	
1,25 ml	25	37,5	50	62,5	75	87,5	100	112,5	125	137,5	
2,50 ml	50	75	100	125	150	175	200	225	250	275	
5,00 ml	100	150	200	250	300	350	400	450	500	550	
10,00 ml	200	300	400	500	600	700	800	900	1000	1100	
12,50 ml	250	375	500	625	750	875	1000	1125	1250	1375	
25,00 ml	500	750	1000	1250	1500	1750	2000	2250	2500	2750	
50,00 ml	1000	1500	2000	2500	3000	3500	4000	4500	5000	5500	
	48	35	23	19	15	13	11	9	8	7	
		Steps									

Contamination-free: Ripette pro and Ritips

The Ripette pro is predestinated for dispensing of contaminated fluids and volatile liquids. Due to the exclusive acceptance of dispenser tips with plunger e.g. Ritips and other compatible tips, the Ripette pro is secure from contamination by ascending aerosols in the tip and guarantees a long-term secure use even with samples of high fluidity.

Assignment of registered trademarks;

Multipette ®: Eppendorf AG, Mobile Phone Step ®: Brand GmbH + Co. KG, Distriman ®: Gilson Inc., bioclean ®: Ritter GmbH



- Due to its low weight the Ripette pro is most-suitable for long dispensing series.
- Ergonomic design which allows an one-hand-operation of volume selection, loading and dispensing of liquids.
- The dosing table shows all dosing possibilities at a glance for choosing the most-suitable tip for your application.
- The Ripette pro is reliable and maintenance-free due to its robust design and a minimum of wear parts.
- By using selected premium materials the Ripette pro shows a high resistance towards a broad range of chemicals.
- Every device will be tested in terms of functionality and precision after production. Only after the requirements of DIN EN ISO 8655 have been fulfilled, a serial number will be issued.



Quasar Instruments, LLC
4835 Centennial Blvd.
Colorado Springs, CO 80919
t. 877-535-8001 | d. 719-535-8000
QuasarInstruments.com