

## A strong start

### MML monocular - compact and inexpensive

Monocular microscopes are ideal for many applications in the laboratory, teaching and production. All models have 45° inclined viewing and 360° rotating optical head. The sturdy metal tripod ensures high stability and the option to choose between coarse and fine adjustment facilitates precise working. The microscopes are equipped with integrated illumination.

Power supply: 230 V (115 V optional).  
A wide range of accessories is available for all models.



MML1200



MML1500

	Optical equipment	Illumination	Special features:	Application
<b>MML1200</b>	10x plane eyepiece Objectives (achromatic) 4x/NA 0.10 10x/NA 0.25 40x/NA 0.65	6 V 10 W brightfield Abbe condenser		Schools Training Simple laboratory applications
<b>MML1300</b>	10x plane eyepiece Objectives (achromatic) 4x/NA 0.10 10x/NA 0.25 40x/NA 0.65	6 V 10 W brightfield Abbe condenser, adjustable	Adjustable illumination	Schools Training Simple laboratory applications
<b>MML1400</b>	10x plane eyepiece Objectives (achromatic) 4x/NA 0.10 10x/NA 0.25 40x/NA 0.65 100x/NA 1.25 oil	6 V 10 W brightfield Abbe condenser, adjustable	Adjustable illumination 4x objective turret	Schools Training Simple laboratory applications
<b>MML1500</b>	10x plane eyepiece Objectives (achromatic) 4x/NA 0.10 10x/NA 0.25 40x/NA 0.67 100x/NA 1.25 oil	6 V 20 W brightfield Abbe condenser, adjustable	Adjustable illumination 4x objective turret XY cross table	Schools Training Simple laboratory applications