



Riplate medio 0,5 ml

Riplate medio 1 ml



Riplate medio 1 ml



Riplate medio 0,5 ml

Riplate® medio

50% weniger Stapelhöhe
Saves 50% of storage capacity

Artikel-Nr.	Artikelbezeichnung	VE	Stück / Umkarton
Article-no.	Article name	PU	Pieces / cover. box
43001 - 0101	Riplate Medio, PP, - 1ml	32	320
43001 - 0102	Riplate Medio, PP, - 0,5ml	40	400

Riplate in den Farben blau, grün, rot, gelb sind auf Anfrage erhältlich.

Riplate deepwell-plates are available in blue, green, yellow and red colour on request.

Riplate medio

- Durch kompakte Bauweise bietet Riplate medio bis zu 50% geringere Stapelhöhe.
- Die erhöhten Ränder der einzelnen Wells lassen die Riplate medio luftdicht verschweißen oder verkleben.
- Die alphanumerische Kennzeichnung ermöglicht die Registrierung und Identifizierung der einzelnen Proben.
- Neu-konzipierte Wells ermöglichen eine gute Probendurchmischung.

Riplate® medio 0,5 ml

- Höhe nur 20,5 mm
- Gesamtvolumen 0,6 ml/well
- Arbeitsvolumen 20 µl – 500 µl

Riplate® medio 1 ml

- Höhe nur 27,0 mm
- Gesamtvolumen 1,2 ml/well
- Arbeitsvolumen 50 µl – 1000 µl

Riplate medio

- *Due to the compact construction of Riplate medio less than 50% stacking size is needed.*
- *Higher well rims allow an airtight closure or sealing of riplate medio.*
- *The alphanumeric marking allows registration and identification of single assays.*
- *New developed conical wells provide a good bioturbation.*

Riplate® medio 0,5 ml


- *Height only 20,5 mm*
- *Total volume 0,6 ml/well*
- *Working volume 20 µl – 500 µl*

Riplate® medio 1 ml

- *Height only 27,0 mm*
- *Total volume 1,2 ml/well*
- *Working volume 50 µl – 1000 µl*



Riplate Plus mit Luftschleusen und
Randanhebung der Wells



Riplate plus



Capmat

Riplate® 1ml

Deepwellplatte im multifunktionelle Standard
Deepwell-plates – standard in multifunction

liquid handling

robotic consumables

clinical products

Artikel - Nr.	Artikelbezeichnung	VE	Stück / Umkarton
Article - no.	Article name	PU	Pieces / cover. box
43001 - 0010	Riplate, PS, - 1 ml	20	200
43001 - 0016	Riplate, PP, - 1 ml	20	200
43001 - 0116	Riplate Plus, PP, - 1 ml	20	200
43001 - 1016	Riplate, PP, - 1 ml, sterile	10	100
43001 - 1116	Riplate Plus, PP, - 1 ml, sterile	10	100
43001 - 0017	Capmat 1ml	20	160
43001 - 1017	Capmat 1 ml, sterile	10	100

Riplate 1 ml

- Erhöhte Wells verringern das Risiko der Kreuzkontamination und verbessern die Sicherheit beim Verschweißen oder Verschließen mit Klebefolie.
- Optimierte Kontur im SBS-Standard und acht Luftschleusen.
- Für Anwendungen im Thermocycler oder zum Auftauen gefrorener Proben im Wasserbad.
- Für Anwendungen mit extremen Temperaturbedingungen wie Tiefgefrieren oder Auftauen in der Mikrowelle.
- Der runde Boden garantiert eine hohe Probenrückgewinnung.
- Durch die chemische Resistenz für sensible Diagnostik-Anwendungen geeignet.

Riplate 1ml

- *Raised wells reduce the risk of cross-contamination and increase the security of hot sealing or sealing with adhesive foils.*
- *Optimized footprint according to the SBS-standard and eight ventilation holes.*
- *Ideal for use in a thermocycler or for defrosting of samples in a hot-water bath.*
- *For freezing or defrosting in a microwave.*
- *The U-shaped well bottom reduces the death-volume.*
- *Due to the high chemical resistance ideal for sensitive diagnostic applications.*